

# Product Information

## Mini Super<sup>HT</sup> Pap Pen Super<sup>HT</sup> Pap Pen

### Catalog Number:

Mini Super<sup>HT</sup> Pap Pen: 22005

Super<sup>HT</sup> Pap Pen: 22006

**Unit Size:** 1 pen

### Size:

Mini Super<sup>HT</sup> Pap Pen: 2.5 mm tip, ~400 applications

Super<sup>HT</sup> Pap Pen: 4 mm tip, ~800 applications

Number of applications varies depending on barrier size

**Storage and handling:** Store at room temperature. The pens contain organic solvents. Avoid skin and eye contact and keep away from open flame. Cap tightly after each use.

### Product Description

Super<sup>HT</sup> Pap Pens are used to create hydrophobic barriers around tissue sections on glass slides to hold staining solutions in place on the slide. This allows the conservation of antibody staining solution or staining of two sections on the same slide with different solutions. The Mini Super<sup>HT</sup> Pap Pen has a fine point that is useful for separating multiple sections on the same slide. Super<sup>HT</sup> Pap Pen barriers are insoluble in aqueous buffers, detergents, alcohol and acetone, but can be removed with xylene. Barriers are stable at temperatures up to 120°C.

### Directions

Before first use, fill the tip with barrier fluid. Press the tip down on a clean glass microscope slide to open the tip valve (Figure 1). Keep pressing down until you see the tip fill completely with green liquid. Stop pressing as soon as the tip is completely filled.

Barriers should be created after deparaffinization and rehydration for paraffin sections. For frozen sections, create barriers when slides are dry, before the first buffer incubation step. To create a barrier, dry the area around the section with a cotton swab if necessary. Draw an unbroken circle of Pap Pen fluid around the section (Figure 2). Do not press the valve down while drawing the barrier. Take care to prevent the Pap Pen fluid from touching the tissue section. Let the barrier dry 15-30 seconds before immersing slides or adding buffer. Dry the area immediately outside the barrier if necessary and add just enough buffer to fill the barrier without overflowing. A square of Parafilm® can be placed on top of the section to spread the buffer evenly and prevent evaporation. For overnight incubations, slides should be placed in a humidified chamber to prevent evaporation.




Figure 1. Filling pen tip before first use. A. Pen tip is packaged dry. B. Press tip down against a microscope slide until tip fills completely with green liquid (C).

### Related Products

Cat.#	Product Name	Unit Size
23001	EverBrite™ Mounting Medium	10 mL
23002	EverBrite™ Mounting Medium with DAPI	10 mL
23003	EverBrite™ Hardset Mounting Medium	10 mL
23004	EverBrite™ Hardset Mounting Medium with DAPI	10 mL
23005	CoverGrip™ Coverslip Sealant	15 mL
40061-T	RedDot™2 Far Red Nuclear Counterstain, 200X in DMSO, Trial Size (15-20 tests)	25 µL
22015	Fixation Buffer	100 mL
22016	Permeabilization Buffer	100 mL
22017	Permeabilization and Blocking Buffer (5X)	100 mL
22010	10% Fish Gelatin Blocking Buffer	100 mL
22011	Fish Gelatin Powder	2 x 50 g
22014	30% Bovine Serum Albumin Solution	100 mL
22002	Tween®-20	50 mL

Please visit [www.biotium.com](http://www.biotium.com) to view our full selection of products featuring bright and photostable fluorescent CF™ dyes, including secondary antibodies and Mix-n-Stain™ antibody labeling kits, and R-PE dye conjugates. Biotium also offers a variety of apoptosis and cell viability assays for flow cytometry analysis, including mitochondrial membrane potential dyes, fluorescent Annexin V conjugates, and NucView™488 Caspase-3 Substrate for live cells.




Figure 2. Creating barriers. A. Draw an unbroken circle around sections with Pap Pen. B. Pipette buffer into the barrier. Add just enough buffer to fill the barrier without overflowing (C).

Parafilm is a registered trademark of Bemis Company. TWEEN is a registered trademark of Uniqema Americas LLC.

Biotium products are high-quality reagents and materials intended for research purposes only. Some products are potentially hazardous chemicals - please read the Material Safety Data Sheet for additional information regarding handling potentially hazardous chemicals. Several of Biotium products and product applications are covered by U.S. and international patents and pending patents. Our products are not available for resale or other commercial uses without a specific agreement from Biotium, Inc. We welcome inquiries about licensing the use of our dyes, trademarks or technologies. Please submit inquiries by e-mail to [btinfo@biotium.com](mailto:btinfo@biotium.com).