

# Product Information

## Mini Cell Scrapers


**Catalog Number:** 22003

**Storage:** Room temperature

**Quantity:** 200 scrapers total, 10 scrapers per sterile packet

### Product Description

Sterile, polyethylene, disposable mini cell scrapers can be used to effectively harvest cells or cell lysates from 24, 48 and 96 well plates. They measure 3/16" (0.5 cm) wide and 2 3/8" (6 cm) long. Mini cell scrapers may be used with the AccuEasy™ Flow Cytometry Kit (catalog number 30069) for lifting highly adherent cells.



### Related Products

Catalog number	Product
30069	AccuEasy™ Flow Cytometry Kit
23006	Flow Cytometry Fixation/Permeabilization Kit
22015	Fixation Buffer
22016	Permeabilization Buffer
22017	Permeabilization and Blocking Buffer
22010	10% Fish Gelatin Blocking Buffer
22011	Fish Gelatin Powder
22014	30% Bovine Serum Albumin Solution
22012	Dry Milk Powder
22002	Tween®-20

Please visit [www.biotium.com](http://www.biotium.com) to view our full selection of products featuring bright and photostable fluorescent CF™ dyes, including secondary antibodies and Mix-n-Stain™ antibody labeling kits. Biotium also offers a variety of apoptosis and cell viability assays for flow cytometry analysis, including mitochondrial membrane potential dyes, fluorescent Annexin V conjugates, and NucView™488 Caspase-3 Substrate for live cells.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.